

Immune Cell Types

Foxp3

Treg Cells

CD3	CD73
CD4	CD152 (CTLA-4)
CD25	GITR
CD39	GARP

ROR- γ

Th17 Cells

CD3	CD196 (CCR6) (h)
CD4	CD161 (h)
CD194 (CCR4) (h)	IL-23R

NK/NKT Cells

CD3	CD57 (h)	CD244
CD11b	CD69	CD314 (NKG2D)
CD11c	CD94	CD319 (h)
CD16	CD122	CD328 (h)
CD25	CD158 (h)	CD335 (NKp46)
CD45R (m)	CD160	CD352 (SLAMF6)
CD49b (m)	CD161 (h)	
CD56 (h)	CD161c (NK1.1) (m)	

Dendritic Cells

CD1a (h)	CD40	CD205 (DEC-205)
CD1b (h)	CD49d (h)	CD207 (Langerin)
CD1c (h)	CD80	CD209 (DC-SIGN)
CD4	CD83	CD273 (B7-DC)
CD8	CD85g (h)	CD289 (TLR9)
CD11b (m)	CD86	CD303 (h)
CD11c	CD123	CD304 (h)
CD14 (h)	CD197 (CCR7)	CD317 (PDCA-1)

Macrophages

CD11a	CD68
CD11b	CD80
CD11c	CD85
CD13	CD86
CD14	CD107a
CD16/CD32	CD115
CD33	CD195
CD64	CD282 (TLR2)
	CD284 (TLR4)

T Cells

CD3	CD44 (m)	CD134
CD4	CD45RA (h)	CD137 (4-1BB)
CD5	CD45RO (h)	CD152
CD7 (h)	CD62L	CD154
CD8	CD69	CD197 (CCR7)
CD25	CD73 (m)	CD272 (BTLA)
CD27	CD101 (h)	CD278 (ICOS)
CD28	CD127	CD279 (PD-1)
CD39	CD130	

Granulocytes

CD11b	CD16
CD11c	CD31
CD13	CD88
CD14	CD114

B Cells

CD2 (m)	CD40	CD93 (m)
CD5	CD45R (m)	CD137 (4-1BB)
CD19	CD69	CD138 (h)
CD20	CD70	CD252
CD21/CD35 (m)	CD79a	CD267 (m)
CD22	CD79b	CD268
CD23	CD80	CD279 (PD-1)
CD25	CD86	

Tonbo Biosciences Research Generics™

The quality of the leading flow cytometry reagents!

Fluorophore-conjugated antibodies, In Vivo Ready™ antibodies, custom and bulk reagents and more!

ANTICUERPOS PARA CITOMETRÍA DE FLUJO.

Los clones que utiliza habitualmente, por un precio mucho más ventajoso.