

Film Products Performance Guide

<u>Imaging subject</u>	<u>Desired Performance</u>	<u>KODAK Film/Products</u>
High energy isotopes ³² P ¹²⁵ I	Maximum sensitivity	BioMax TranScreen HE/BioMax MS film
	High resolution and high sensitivity	BioMax TranScreen LE/BioMax MS film
	Maximum resolution For microscopic examinations	BioMax MR film Emulsion types NTB-2 and NTB-3
Low and medium energy isotopes ³ H, ¹⁴ C, ³⁵ S, ⁴⁵ Ca, ³³ P and ⁵⁹ Fe	Maximum sensitivity	BioMax TranScreen LE/BioMax MS film
	Maximum resolution	BioMax MR film
Fluorographic reagents	Maximum sensitivity	BioMax TranScreen LE/BioMax MS film
	Maximum resolution	BioMax TranScreen LE/BioMax MS film
Chemiluminescent markers	Maximum sensitivity	BioMax Light film
	Maximum resolution	
Electron microscopic images	Maximum image quality	EM film 4489
	High image quality and signal-to-noise ratio	EM film SO-163